

Product information

SENDOLINE AB

H-File/ Hedström File

The cutting ability of the H-File is very good.

The cylindrical core design making the Sendoline H-File one of the most flexible H-Files available today.

It is designed with “one” cutting edge with a 90° tip and a 90° angle. Its unique spiral design flutes successively increase in depth and length from tip to handle, enhancing its efficiency while maintaining safety.

The H-file should be used with a filing up and down movement.

The H-File is also available in nickel titanium named NiTi-H File

ISO/Sizes: 06-140

Length: 21 mm, 25 mm, 28 mm, 31 mm

Package: Individual or assorted (in sizes 15-40 or 45-80) in blister. All packages come with 6 pcs/blister

Material: Operating part: Stainless steel. Handle: Plastic. Stopper: Silicone.

