



## TePe EasyPick™

<b>Product</b>	Plastic core with cleaning surfaces of silicone.
<b>Autoclave</b>	Program: Standard 134°C. Process duration: 55 min, sterilization duration: 7 min, sterilization pressure: 3 bar / 300 kPa.
<b>Pack</b>	12 pcs/clam pack, 36 pcs/clam pack, 10 clam packs/cartoon. The clam pack (blister) is made of Polyethylene terphthalate (PET). The backing card is made of cardboard and the folder of paper. The carton is made of chlorinefree recyclable cardboard. The pocket case is made of polypropylene (PP).
<b>Recycling</b>	Product in the usual household waste. The blister pack should be separated and recycled. The plastic pack and pocket case in container for plastic waste and the backing card in container for recycling of cardboard. The cartons in container for recycling of corrugated cardboard.
<b>Product features</b>	<ul style="list-style-type: none"><li>- Easy cleaning between your teeth.</li><li>- 2 sizes to fit both narrow and wider interdental spaces.</li><li>- Ergonomic handle for easy use, non-slip grip.</li><li>- Efficient silicone cleaning surfaces.</li><li>- Durable.</li></ul>
<b>Developed in collaboration with dental experts</b>	Developed in collaboration with dental experts.